MassageGun Advanced

Art.nr. 373 815
Model TO-18

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Important information

Danger: To reduce the risk of electrical shock, fire or personal injuries this product must be used in accordance with following instructions.

The use of this product is for private use only. Any warranties will be invalid if it is used in commercial or multi-user environments such as sports teams or companies.

- 1) This product can be used for maximum one hour at the time. Then let the product rest for 30 minutes before using it again.
- 2) The product can only be used on and by adults. Adults with disabilities must be supervised by another adult.
- 3) The product should only be used on a dry and clean part of the body.
- 4) When massaging, press lightly over the muscles for about 60 seconds per area.
- 5) Stop all use if you experience some sort of pain or discomfort. Avoid massaging the head or bones.
- 6) Bruising can occur after use regardless of massage level and pressure.
- 7) To avoid squeezing injuries, keep fingers, hair or other small bodyparts away from the massage head holder.
- 8) Keep the ventilation holes free of dust/debris and do not place objects into the ventilation holes.
- 9) The product should not be used outdoors.
- 10) Keep the product away from all types of liquids and heat sources.
- 11) Always inspect the product carefully before use. The product should not be used if it has fallen to the ground, has a damaged cord, a damaged plug or has been in water.
- 12) Do not remove any screws or attempt to disassemble the product.
- 13) The product should only be charged with the included charger.
- 14) Never leave the product switched on or charge unattended.
- 15) Do not use the product without the approval of a doctor.

Reservations are made regarding printing errors in the manual.

Use of massage heads

To attach or remove the massage head, gently push in/pull out the massagehead while twisting.

Round massage head

The round massage head is the most versatile of the 5 massage heads. The round shape provides a good blood flow and a deep massage of the muscles. This makes it excellent for muscle relaxation both before and after exercise. The massage head is recommended to massage large muscles like the abdominal and side muscles abdomen, chest, neck, broad back muscles, buttocks and thighs.

Flat massage head

The flat shape gives a lighter massage of the muscles and is well suited for massage of the straight abdominal, abdominal side muscles, chest, broad back muscles, hips, buttocks, thighs and legs.

U-shaped massage head

The U-shaped massage head is designed to avoid the spine and cervical vertebrae. This makes the shape suitable for the neck, shoulders, and paraspinal muscles. Avoid direct contact with the spine or cervical vertebrae as this can lead to physical injury.

Conical massage head

The conical shape is excellent for a deep massage of the soles of the feet, palms, rhomboids and the broad back muscles. The conical shape can also be used for a deeper massage of muscle tension in thigh and glutes when the round or flat massage head does not reach all the way in.

Thumbshaped massage head

This shape is designed to give a gentle massage. The thumbshaped massage head is recommended to use in palms, underarms and soles of the feet. The thumb shape gives an relaxing feeling and increase the blood flow.

Use, charging and maintenance

Charging

1. To charge the product, attach the included charger to the battery charging port.
2. The battery can be charged at any time and at any battery level.
3. The battery is fully charged after 3,5 hours.
4. Average use time of a fully charged battery is 3+ hours depending on the massage level and pressure used.

Use

1. Press and hold the on/off button in front for 3 seconds to turn on the product.
2. Select the massage head you want to use and press the on/off button in front or the start button on the back to start the massage.
3. There is six different massage levels you can choose from by pushing the start button on the back or on/off button in front multiple times.
4. Massage the desired body part at the desired massage level and apply pressure as needed.
5. When the light at the front turns red it means that there is only 20% battery left. Then we recommend that you charge the product.
6. To turn off the product, push and hold the on/off button in front for 3 seconds.

Maintenance

1. Turn off the product before cleaning.
2. Use a slightly damp cloth to clean the product and wipe with a soft, dry cloth after.

WARRANTY REGULATIONS

The warranty regulations apply for products that are used for home exercise. For the corporate market, there are specially adapted products and separate warranty regulations apply.

Exercise equipment consist mainly of three component categories; ***Wearing parts***, ***Other replaceable labour parts***, and ***Frame / Structure***.

By **Wearing parts** is implied parts that will eventually be worn down through normal use, and thereby need to be replaced regularly. Included in this category will come breaking pads, wires, pulleys, adjustment screws, ball bearings, rubber bands, springs etc.

Normally, there is no statutory right to bring forth claims for such parts, but we offer a 1 year warranty for these parts.

By **Other replaceable labour parts** is implied all other replaceable parts on the machine or device. Examples of such parts could be seat cushions, foam padding, pedals, exercise computers and other electronics. For the consumer market, a 2 year restricted warranty applies for these parts, provided that they have not been subjected to negligence or error in handling or storing, wrong assembly, unsuitable use or unauthorized or improper repairs.

Excluded from this restriction are treadmill motors.

For all treadmill motors, there is a 5 year warranty, provided that they have not been subjected to negligence in maintenance, or unauthorized or improper repairs.

By **Frame/structure** is implied all the parts of the load-bearing steel frame, and its welding seems and joints.

For these parts, the warranty time for a minimum of 5 years is applied. Longer warranties and/or life time warranties will specified in the product information for each individual machine or device.

IMPORTANT REGARDING SERVICE

In the event of problems of any kind, please contact MyIna Service. We would like you to contact us before contacting the store so we can offer you the best possible help.

Visit our website www.mylinasport.com - here you will find information about the products, user manuals, a contact form to get in touch with us and possibility to order spare parts. By filling out the contact form you give us the information we need to help you as effectively as possible.

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Mylna Sport AB